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1. Angelo, Troedhan et al. "Biomechanical Stability of Dental Implants in Augmented Maxillary Sites: Results of a Randomized Clinical Study with Four Different Biomaterials and PRF and a Biological View on Guided Bone Regeneration." BioMed Research International 2015 (2015): 850340. PMC. Web. 25 May 2015.
2. Vivek Gupta, Vivek K. Bains, G. P. Singh, Ashish Mathur, Rhythm Bains, Regenerative Potential of Platelet Rich Fibrin In Dentistry: Literature Review Asian Journal of Oral Health & Allied Sciences -Volume 1, Issue 1, Jan-Mar 2011
3. Marco del corso, Michael Toffler, David M. Dohan, Use of an Autologous Leukocyte and Platelet-Rich Fibrin (L-PRF) Membrane in Post-Avulsion Sites: An Overview of Choukroun's PRF. The journal of implant & Advanced Clinical Dentistry. 27-35. Vol.1, No.9. December/January 2010
4. Tatullo M, Marrelli M, Cassetta M, et al. Platelet Rich Fibrin (P.R.F.) in Reconstructive Surgery of Atrophied Maxillary Bones: Clinical and Histological Evaluations. International Journal of Medical Sciences. 2012;9(10):872-880. doi:10.7150/ijms.5119.
5. Sánchez Garcés MA, Pié Sánchez J, Harutunian K, Forri A, González Martínez R, Real Osuna J, Fierro Garibay C, Gay Escoda C. Revisión bibliográfica de Implantología Bucofacial del año 2008.
6. Naik B, Karunakar P, Jayadev M, Marshal VR. Role of Platelet rich fibrin in wound healing: A critical review. Journal of Conservative Dentistry. JCD. 2013;16(4):284-293.

ANKYLOS® | DENTSPLY IMPLANTS